



CE 1275	Nylon -Polyamide 6/6 Spatula / lancet tip						Made in Germany	
	USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE		ART. CODE
SPATULA	11/0	15	3/8	4.2	100	S	A007PA	
		30	1/2	5.0	150		A090PA	
		30	1/2	5.6	150		A088PA	
	10/0	15	3/8	4.2	100	S	A008PA	
		30		4.2	100		A073PA	
		15		4.8	150	S	A052PA	
		30		4.8	150	S	A009PA	
		45		4.8	150	S	A058PA	
		30		4.8	150		A010PA	
		38		4.8	150		A047PA	
		30		4.8	200		A067PA	
		45		4.8	200		A055PA	
		15		5.6	150	S	A075PA	
		30		5.6	150	S	A061PA	
		30		5.6	150		A051PA	
		15		6.0	100		A037PA	
		30		6.0	100		A039PA	
		15		6.0	150	S	A092PA	
		30		6.0	150		A062PA	
		38		6.0	150		A066PA	
		30		6.0	200		A102PA	
		4		6.2	150	S	A024PA	
		15		6.2	150	S	A043PA	
		20		6.2	150	S	A013PA	
		30		6.2	150	S	A014PA	
		38	3/8	6.2	150	S	A069PA	
		10		6.2	150		A087PA	
		15		6.2	150		A085PA	
		20		6.2	150		A034PA	
		30		6.2	150		A016PA	
		38		6.2	150		A036PA	
		45		6.2	150		A045PA	
		75		6.2	150		A048PA	
		15		6.4	150	S	A091PA	
		30		6.4	150		A080PA	
		45		6.4	150		A065PA	
		15		6.4	200	S	A076PA	
		30		6.4	200	S	A025PA	
		38		6.4	200	S	A020PA	
		45		6.4	200	S	A074PA	
		30		6.4	200		A026PA	
		38		6.4	200		A023PA	
		45		6.4	200		A041PA	
		75		6.4	200		A053PA	
		38		6.4	220		A057PA	
		38		6.6	200		A060PA	
		30	3/8	7.0	200		A070PA	
		15	1/2	5.2	150		A098PA	
		15		5.6	150	S	A040PA	
		20		5.6	150	S	A001PA	
		30		5.6	150	S	A002PA	
		30		5.6	150		A004PA	
		30		6.0	150		A071PA	
		13		6.2	100	S	A054PA	
		30		6.2	150		A101PA	
		23		6.4	200	S	A081PA	
		45		6.4	200		A068PA	
		30		7.1	150		A096PA	
		30	1/2	8.0	200		A077PA	
		15	50-90°	4.8	150	S	A029PA	BICURVED
		20		4.8	150		A031PA	
		30		4.8	150		A032PA	
		30		5.0	150		A089PA	
		30		6.0	200		A102PA	
		38	50-90°	6.4	200		A104PA	BICURVED
		30	Straight	16.0	150	S	A097PA	


CE 1275	Nylon -Polyamide 6/6 Spatula / lancet tip						Made in Germany
	USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE	

SPATULA



NYLON

9/0	30	1/4	6.4	200		A027PA	
	45		6.5	200		A083PA	
	45		6.0	200		A084PA	
	30		6.4	200		A028PA	
	45		8.0	300		A063PA	
	45		8.5	300		A064PA	
	45		8.5	450		A099PA	
	75	1/4	8.5	450	S	A100PA	
	15	3/8	4.8	150	S	A030PA	
	30		4.8	150		A011PA	
	15		6.0	100		A038PA	
	30		6.0	150		A072PA	
	20		6.2	150	S	A015PA	
	30		6.2	150	S	A079PA	
30	6.2		150		A017PA		
30	6.4		150		A093PA		
30	6.4		200	S	A018PA		
38	6.4		200	S	A044PA		
30	6.4		200		A021PA		
38	6.4		200		A042PA		
45	6.4		200		A049PA		
75	6.4		200		A050PA		
75	3/8	8.0	200		A086PA		
20	1/2	4.8	150		A094PA		
20		5.6	150	S	A003PA		
20		5.6	150	S	A003PA		
30		5.6	150		A005PA		
30		7.1	150		A056PA		
30		1/2	8.0	200		A078PA	
20	50-90°	4.8	150		A033PA	BICURVED	
8/0	30	3/8	4.8	150		A012PA	
	30		6.4	200	S	A019PA	
	38		6.4	200	S	A035PA	
	45		6.4	200	S	A059PA	
	30		6.4	200		A022PA	
	38		6.4	200		A046PA	
	30		6.4	200		A022PA	
	38	3/8	6.4	200		A046PA	
30	1/2	5.6	150		A006PA		
7/0	30	3/8	6.4	200		A095PA	
	23	1/2	10.0	230		A082PA	
6/0	45	1/4	8.0	300		A063PA	
5/0	45	1/4	8.5	350		A064PA	
	45	3/8	10.0	350		A103PA	
4/0	45	1/4	8.5	450		A099PA	
	75	1/4	8.5	450	S	A100PA	

CE 1275	Braided Silk -Spatula / Lancet tip						Made in Germany		
	USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE		ART. CODE	
 <p>Braided Silk</p>	10/0	15	3/8	4.2	100	S	A008SB		
		30		4.8	150	S	A009SB		
		30		4.8	150		A011SB		
		15		5.6	150	S	A049SB		
		30		6.0	150		A051SB		
		30		6.2	150	S	A052SB		
		45		6.2	150		A044SB		
		45		6.4	200		A045SB		
		30		6.4	200	S	A053SB		
		30		6.4	200		A054SB		
		20		1/2	5.6	150	S	A001SB	
		30		1/2	5.6	150		A004SB	
		9/0		45	1/4	6.4	200		A021SB
				30	3/8	4.8	150		A012SB
30	6.0		150		A050SB				
30	6.4		200	S	A017SB				
30	6.4		200		A019SB				
38	6.4		200		A039SB				
45	3/8		6.4	200		A041SB			
20	1/2		5.6	150	S	A002SB			
30	5.6		150		A005SB				
30	1/2		7.1	150		A038SB			
8/0	45	1/4	6.4	200		A022SB			
	30	3/8	4.8	150		A013SB			
	30	6.0	150		A047SB				
	30	6.0	150		A065SS				
	45	6.0	150	S	A067SS				
	30	6.2	150	S	A014SB				
	30	6.2	150		A015SB				
	38	6.2	150		A036SB				
	45	6.2	150	S	A010SB				
	45	6.2	150		A016SB				
	30	3/8	6.4	200	S	A018SB			
	30	6.4	200	S	A061SB				
	30	6.4	200		A020SB				
	38	6.4	200	S	A029SB				
	38	6.4	200		A032SB				
	45	6.4	200	S	A027SB				
	45	6.4	200		A026SB				
	45	6.4	200		A057SS				
	75	6.4	200		A042SB				
	75	8.0	200		A034SB				
	75	3/8	9.7	230		A031SB			
	20	1/2	5.6	150	S	A003SB			
	30	5.6	150		A006SB				
	30	8.0	220		A066SS				
	45	1/2	8.0	250		A060SS			
7/0	45	3/8	6.4	200		A046SS			
	75	6.4	200		A030SB				
	30	8.0	250		A059SB				
	75	3/8	8.0	250		A058SS			
	45	1/2	8.0	250		A033SS			
6/0	45	1/4	6.4	350		A035SS			
	45	8.0	340		A064SS				
	45	1/4	8.5	350		A024SS			
	45	3/8	6.4	200		A048SS			
	45	8.0	300		A040SS				
	75	8.0	300		A055SS				
	75	8.0	340		A043SB				
	45	3/8	8.5	350		A028SS			
	45	1/2	6.0	300		A037SS			
	45	1/2	8.5	350	S	A007SS			
5/0	20	1/4	8.0	460	S	A069SS			
	45	1/4	8.5	350		A025SS			
	45	3/8	8.0	350	S	A056SS			
	20	3/8	8.0	460	S	A068SS			
4/0	30	1/4	6.0	350		A023SB			
	20	1/4	8.0	460	S	A063SS			
	20	3/8	8.0	460	S	A062SS			

CE 1275	PGA -Polyglycolic Acid Spatula / lancet tip Violet /Undyed						Made in Germany
	USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE	



PGA

USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE	ART. CODE	Made in Germany
10/0	15	Bicurved	4.8	150	S	A062PGA	50-90°
	10	3/8	4.8	150	S	A035PGA	
	30		6.0	150		A061PGA	
	10		6.2	150	S	A002PGA	
	30		6.2	150		A083PGA	
	45	3/8	6.2	150		A045PGA	
	30	1/2	5.5	220	S	A058PGA	
	30		5.6	150		A079PGA	
	30		6.0	150	S	A070PGA	
	15		6.2	150	S	A042PGA	
	30	1/2	6.2	150	S	A066PGA	
	9/0	30	3/8	6.4	200		
30			6.4	200	S	A044PGA	
30			6.4	200		A046PGA	
30		3/8	6.4	200	S	A063PGA	
30		1/2	5.5	220	S	A059PGA	
30			6.0	150	S	A071PGA	
30			6.0	150		A074PGA	
30			6.2	150	S	A067PGA	
30		1/2	6.2	150		A068PGA	
8/0		70	3/8	5.6	150		A004PGA
	15	3/8	6.0	200	S	A049PGA	
	30	Bicurved	6.0	200	S	A033PGA	
	45	3/8	6.0	200		A027PGA	
	30		6.2	150		A043PGA	
	30		6.4	200	S	A051PGA	
	30		6.4	200		A003PGA	
	45		6.4	200		A030PGA	
	45		6.4	200		A082PGA	
	30		8.0	220		A016PGA	
	45	3/8	8.0	200		A050PGA	
	30	1/2	8.0	220	S	A012PGA	
	30	1/2	8.0	250		A060PGA	
	7/0	45	3/8	6.0	220		A014PGA
30			6.4	200	S	A065PGA	
30			6.4	220		A011PGA	
30			6.4	200		A053PGA	
45			6.4	200	S	A009PGA	
45			6.4	200		A034PGA	
30			8.0	220		A015PGA	
45			8.0	230		A005PGA	
45			9.6	230		A080PGA	
45		3/8	10.0	250	S	A039PGA	
30		1/2	6.4	200		A064PGA	
45		1/2	8.0	250		A038PGA	

CE 1275	PGA -Polyglycolic Acid Violet /Undyed						Made in Germany
	USP	LENGTH (cm)	CURVE	LENGTH (mm)	DIAM. (micron)	ARMED S- SINGLE	



PGA

6/0

45	1/4	6.6	250		A054PGA
45		8.0	300		A001PGA
75		8.0	300		A032PGA
45		8.0	330		A072PGA
45		8.0	340		A020PGA
45		8.0	350		A006PGA
30		8.0	350		A018PGA
45		8.0	350		A036PGA
45		8.5	350		A057PGA
45		9.0	350		A040PGA
45		9.5	350		A077PGA
45	1/4	10.0	340		A048PGA
45	3/8	6.4	200		A052PGA
45		8.0	300		A024PGA
75		8.0	300		A031PGA
45		10.0	340		A025PGA
75	3/8	10.0	340		A026PGA
30	1/2	6.4	220		A069PGA
30		6.4	240		A075PGA
45	1/2	8.0	300		A041PGA

5/0

45	1/4	8.0	300		A076PGA
45		8.0	340		A019PGA
45		8.0	340		A055PGA
75		8.0	350	S	A010PGA
30		8.0	350		A017PGA
45		8.0	380		A073PGA
45		8.5	350		A007PGA
45		8.5	350		A037PGA
45	1/4	10.0	340		A047PGA
45	3/8	8.0	340		A021PGA
75		8.0	340		A028PGA
30		8.0	350	S	A008PGA
45		10.0	300	S	A013PGA
45		10.0	340		A022PGA
75	3/8	10.0	340		A023PGA
45	1/2	8.0	350		A081PGA

4/0

45	1/4	8.0	440		A056PGA
45	3/8	8.0	450		A078PGA